

ELECTRIC HEATING INCUBATOR



HPX series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

ELECTRIC HEATING INCUBATOR >>

HPX series

Product Features >

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of mirror stainless steel, the shelves can be adjusted optionally;
- Adopting PMMA I operating system, temperature control, timing, and over temperature alarm;
- Imported motor and fan, intelligent temperature controller with timing function;
- Observation window on the door, convenient for users to observe;
- Nano material magnetic door seal strips and heat preservation material;
- Power supply recovery function, data will not be lost when power failure and system halted.



Shelf



Temperature Sensor



Control Panel



Model	HPX-9052MBE	HPX-9082MBE	HPX-9162MBE	HPX-9272MBE
Volume [L]	40	76	142	234
Display	LED			
Temp. Range [°C]	RT+5~65			
Temp. Resolution [°C]	0.1			
Temp. Fluctuation at 37°C [°C]	±0.5			
Temp. Uniformity at 37°C [°C]	±1.5			
Power Supply [V, Hz]	220±10%, 50±2%			
Power [W]	300	450	450	750
Timer Range(min/h)	0~999			
Internal Dimensions W×D×H [mm]	350×330×350	450×380×450	550×470×550	600×520×750
Housing Dimensions W×D×H [mm]	640×556×505	740×606×605	840×696×705	890×746×905
Number of Shelves (std./max.) [pc]	2/4	2/5	3/7	3/10
Net Weight /Gross Weight [kg]	37/47	50/61	68/91	88/111
Package Size [mm]	740×696×645	840×746×745	960×856×855	1010×906×1055