

ILLUMINATION INCUBATOR



BXL Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

ILLUMINATION INCUBATOR >>

BXL Series

Product Features >

- BRIGHT I control system, automatically compensate parameters based on the environment;
- 4.3 inch LCD touch screen;
- Support for viewing alarm records, data, and curves;
- Supports storage of 100000 records and USB data export, multi-level password management;
- Variable frequency cooling, adjustable height partition type light source;
- Full spectrum LED light, each layer of lighting panel can be individually adjusted.



Doorknob



LED Light



Shelf

Model	BXL-250	BXL-400	BXL-80	BXL-1000
Volume [L]	250	400	800	1000
Temp. Range [°C]	0~70	0~70	0~70	0~70
Temp. Resolution [°C]	0.1	0.1	0.1	0.1
Temp. Fluctuation at 25°C [°C]	±0.5	±0.5	±0.5	±0.5
Temp. Uniformity at 25°C [°C]	±1	±1	±1.5	±1.5
Illuminance [Lux]	0-40000	0-40000	0-40000	0-40000
Illuminance [Lux]	0-40000	0-40000	0-40000	0-40000
Power [W]	1450	1750	3100	3350
Growth Height [cm]	45	45	45	45
Total Growth Area [m ²]	0.46	0.77	1.43	1.83
Internal Dimensions W×D×H [mm]	510×450×1090	600×640×1050	1220×585×1123	1363×670×1123
Housing Dimensions W×D×H [mm]	650×740×1726	745×930×1695	1503×910×1820	1505×1006×1811
Number of Shelves (std./max.) [pcs]	2/4	2/4	2/4	2/4
Number of Layer Illumination Shelves (std./max.) [pcs]	2/4	2/4	2/4	2/4
Net Weight [kg]	143	169	330	340
Gross Weight [kg]	173	201	386	390
Package Size [mm]	790*930*1930	890*1110*1900	1640*1090*2040	1640*1180*2040